EN ELECTRODES FOR DEFIBRILLATION, PACING, MONITORING, CARDIOVERSION

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93465 10078494_a 27.06.2023

EN

ELECTRODES FOR DEFIBRILLATION, PACING, MONITORING, CARDIOVERSION (for adults and pediatrics/infant care)
Please read this document and keep it safe. Ensure that any individual using this product is familiar with and understands the information contained in this document.

Intended Use / Intended Purpose
Defibrillation electrodes / multifunction electrodes are provided as an accessory for automatic external defibrillators (AED), semi-automatic defibrillators or manual defibrillators.
Defibrillation electrodes / multifunction electrodes are applied to the patient's upper body. Defibrillation electrodes / multifunction electrodes are non-sterile, intended only for single use on intact (uninjured) skin and are used to transfer energy or ECG signals between the device and the patient.

To improve readability, the defibrillation electrodes / multifunction electrodes are hereinafter referred to as "electrodes" in this document.

WARNING
- Not following these instructions or any other misuse or misapplication of any electrode may result in patient burns or ineffective therapy. (Note: Reddening of skin is normal.)
- Do not operate the electrodes in the immediate presence of flammable gases (including concentrated oxygen) and solvents to avoid possible explosion or fire hazard. Avoid contact between the product and solvent-based liquids.
- Use these electrodes only with compatible defibrillators as listed on the packaging. Use of this product with any device not listed may result in potentially hazardous connection.
- Strangulation Hazard – Route cables away from patient's throat.
- Non-invasive pacing may cause skin irritation or burns. Use with care and inspect skin every 30 minutes during continuous pacing. Replace electrodes after specified time of continuous pacing (see below "For adult patients" and "For pediatric patients/infants"). Extensive use of the electrodes may lead to air bubbles. To ensure optimal contact with the skin, smooth out bubbles by hand.
- Use only electrodes that are suitable and appropriately labeled for X-ray imaging. Unlike radiotransparent electrodes, radiolucent electrodes may lead to shadows on X-ray images. The cables and cable connections block X-rays and therefore appear as white patches on X-ray images. The physician needs to check in advance if the radiolucency of the radiolucent electrodes is sufficient for application.

PRODUCT LIMITATIONS

For adult patients:
- For use only on adult patients.
- Replace electrodes after:
 - 50 defibrillation shocks
 - 1 hour of continuous pacing using
 140 mA / 120 bpM / 20 ms pulse duration
 80 mA / 100 bpM / 40 ms pulse duration
 - 8 hours of continuous pacing using
 70 mA / 60 bpM / 20 ms pulse duration
 50 mA / 60 bpM / 40 ms pulse duration
 - 24 hours on skin

For pediatric patients/infants:
- For use only on pediatric patients/infants
- For use only with manual defibrillators or in manual mode. Do not use with automatic external defibrillators (AEDs), unless the AED model is explicitly listed on the packaging as compatible.
- Replace electrodes after:
 - 25 defibrillation shocks
 - 1 hour of continuous pacing using
 70 mA / 140 bpM / 20 ms pulse duration
 70 mA / 140 bpM / 40 ms pulse duration
 - 8 hours on skin

COMPATIBILITY
Use these electrodes only with compatible defibrillator devices as listed on the packaging. Use of this product with any device not listed may result in potentially hazardous connection.

CAUTION
- Do not use once product has expired.
- Do not crush, bend or fold electrodes or store them under heavy objects.
- Do not open packaging until ready for use.
- Use electrodes within 1 day once packaging has been opened.
- Inspect packaging, electrode, cable and connector prior to use. Do not use if compromised.
- Do not use if gel is dry.
- Do not use additional gel on electrodes.
- Ensure patient's chest is clean and dry.
- Do not overlap the electrodes.
- Keep applied defibrillation electrodes clear of any other electrodes or metal parts in contact with patient. Avoid placing electrodes above an internal pacemaker.
- Do not reposition electrodes.
- Use separate ECG electrodes and cable when performing non-invasive pacing.
- Do not touch patient during defibrillation.
- Do not discharge hand-held paddles through these electrodes.
- Do not re-use the electrodes.
 If re-used, the electrical properties may be insufficient, which could lead to a patient injury/hazard.
- These electrodes are single-patient use. Reuse results in a risk of cross-infection from one patient to another.
- Always keep a second pack of electrodes with the defibrillator.
- For MRT (magnetic resonance tomography) applications, only use appropriately labeled electrodes.

INSTRUCTIONS FOR USE
The electrodes must only be used by appropriately qualified medical personnel, who have been trained in how to use them based on these instructions for use. This does not apply to users of AEDs (automatic external defibrillators), for which instruction is given via voice prompts from the defibrillator itself and/or visual instructions on the defibrillator's display.

NOTE
- The operating instructions of the device in use must also be observed.
- Use local protocol, American Heart Association or European Resuscitation Council guidelines to select appropriate defibrillation energy level.

APPLICATION OF DEFIBRILLATION ELECTRODES / MULTIFUNCTION ELECTRODES
1. Select application sites on the patient's exposed upper body according to the defibrillator manufacturer's operating instructions. There are diagrams to show the two standard options; for adult patients and for pediatric patients.
 - (3a) Anterior-Anterior electrode placement (for defibrillation).
 - (3b) Anterior-Posterior electrode placement (for defibrillation and pacing).
2. Ensure the application site is clean and dry (1/1b). Shave if required (1a).
3. Tear open the packaging (2a), remove electrodes from packaging and then remove the protective cover from the electrodes (2b/2c).
4. Apply to patient in accordance with protocol. Press down firmly (4) to ensure good contact of the entire adhesive surface with the skin.
5. If the electrode is not pre-connected, connect the electrodes to the defibrillator (5).
6. Proceed with defibrillation/pacing according to defibrillator manufacturer's operating instructions.
 Make sure that the electrodes are not used for longer than is stated on the packaging (6).
7. Disconnect electrode connector from defibrillator.
8. After use, carefully remove the electrode using one hand whilst supporting the skin beneath it with the other. Lift the electrode at the end where the electrode cable is attached and slowly

peel it off. Ripping off the electrode or peeling it off quickly may damage the skin. Exercise particular caution if the skin is exceptionally sensitive, especially in children but also in elderly patients, diabetics, or patients on a prolonged course of certain medications known to cause dermatological side effects.

9. Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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