corPatch electrodes

EN ELECTRODES FOR DEFIBRILLATION, PACING, MONITORING, CARDIOVERSION

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93502 10080359_a 14.09.2023

EN

ELECTRODES FOR DEFIBRILLATION, PACING, MONITORING AND CARDIOVERSION (for adults and pediatrics)

Please read this document and keep it safe. Ensure that any individual using this product is familiar with and understands the information contained in this document.

Intended Use

Defibrillation electrodes / multifunction electrodes are provided as an accessory for automatic external defibrillators (AED), semi-automatic defibrillators or manual defibrillators.

Defibrillation electrodes / multifunction electrodes are applied to the patient's upper body. Defibrillation electrodes / multifunction electrodes are non-sterile, intended only for single use on intact (uninjured) skin and are used to transfer energy or ECG signals between the device and the patient.

To improve readability, the defibrillation electrodes / multifunction electrodes are hereinafter referred to as "electrodes" in this document.

WARNING

- Not following these instructions or any other misuse or misapplication of any electrode may result in patient burns or ineffective therapy. (Note: Reddening of skin is normal.)
- Do not operate the electrodes in the immediate presence of flammable gases (including concentrated oxygen) and solvents to avoid possible explosion or fire hazard.
 Avoid contact between the product and solvent-based liquids. Electrodes REF 04324.3, REF 05120.3, REF 05120.4 and REF 05120.5 are approved for operation in multiplace hyperbaric chambers up to 3 barg and an oxygen concentration of <23%.
- Only use the electrodes with **corpuls** defibrillators. Use of this product with any other device may result in a potentially hazardous connection.
- Strangulation hazard – route cables away from patient's throat.
- Non-invasive pacing may cause skin irritation or burns. Use with care and inspect skin every 30 minutes during continuous pacing. Replace electrodes after the specified period of continuous pacing (see below "For adult patients" and "For pediatric patients/infants"). Extensive use of the electrodes may lead to air bubbles. To ensure optimal contact with the skin, smooth out bubbles by hand.

PRODUCT LIMITATIONS

For adults and pediatric patients
(REF 04324.3, REF 05120.4, REF 05120.5):

- Only for use on adults and pediatric patients. Replace the electrodes after:
- Replace electrodes after:
 - ► 50 defibrillation shocks
 - ► 1 hour of continuous pacing using
 140 mA / 120 bpM / 22 ms pulse duration
 80 mA / 100 bpM / 40 ms pulse duration
 - ► 8 hours of continuous pacing using
 70 mA / 60 bpM / 22 ms pulse duration
 50 mA / 60 bpM / 40 ms pulse duration
 - ► 24 hours on skin

For pediatric patients/infants (REF 05120.3):

- Use only on pediatric patients/infants.
- Replace electrodes after:
 - ► 25 defibrillation shocks
 - ► 1 hour of continuous pacing using
 70 mA / 140 bpM / 22 ms pulse duration
 70 mA / 140 bpM / 40 ms pulse duration
 - ► 8 hours on skin

COMPATIBILITY

Use these electrodes only with compatible defibrillator devices as listed on the packaging. Use of this product with any device not listed may result in potentially hazardous connection.

- Use REF 04324.3, REF 05120.3, REF 05120.4 and REF 05120.5 with adapter cable REF 02125.3CA when using **corpuls** 08/16 defibrillator.
- Use the adapter cable specified on the packaging.
- Corpuls AED: Only compatible with REF 05120.3 and REF 05120.4

CAUTION

- Do not use once product has expired.
- Do not crush, bend or fold electrodes or store them under heavy objects.
- Do not open packaging until ready for use.
- Use electrodes within 1 day once packaging has been opened.
- Inspect packaging, electrode, cable and connector prior to use. Do not use if compromised.
- Do not use if gel is dry.
- Do not use additional gel on electrodes.
- Ensure patient's chest is clean and dry.
- Do not overlap the electrodes.
- Keep applied defibrillation electrodes clear of any other electrodes or metal parts in contact with patient. Avoid placing electrodes above an internal pacemaker.
- Inspect pads after repositioning electrodes for performing non-invasive pacing. Do not use if electrodes are damaged.
- Use separate ECG electrodes and cable when performing non-invasive pacing.
- Do not touch patient during defibrillation.
- Do not discharge hand-held paddles through these electrodes.
- Do not re-use the electrodes.
 If re-used, the electrical properties may be insufficient, which could lead to a patient injury/hazard.
- These electrodes are single-patient use. Reuse results in a risk of cross-infection from one patient to another.
- Always keep a second pack of electrodes with the defibrillator.
- For MRT (magnetic resonance tomography) applications, only use appropriately labeled electrodes.

INSTRUCTIONS FOR USE

The electrodes must only be used by appropriately qualified medical personnel, who have been trained in how to use them based on these instructions for use. This does not apply to users of AEDs (automatic external defibrillators), for which instruction is given via voice prompts from the defibrillator itself and/or visual instructions on the defibrillator's display.

NOTE

- The operating instructions of the device in use must also be observed.
- Use local protocol, American Heart Association or European Resuscitation Council guidelines to select appropriate defibrillation energy level.

APPLICATION OF DEFIBRILLATION ELECTRODES / MULTIFUNCTION ELECTRODES

1. Select application sites on the patient's exposed upper body according to the defibrillator manufacturer's operating instructions. The illustrations show the two standard options; for adult patients and for pediatric patients.
 - (3a) Anterior-Anterior electrode placement (for defibrillation)
 - (3b) Anterior-Posterior electrode placement (for defibrillation and pacing)
2. Ensure the application site is clean and dry (1/1b). Shave if required (1a).
3. Tear open pouch or foil (2a), remove electrodes from packaging and then remove the protective foil from the electrodes (2b).
4. Apply to patient in accordance with protocol (3a/3b). Press down firmly (4) to ensure good contact of the entire adhesive surface with the skin.
5. Connect the electrodes to the defibrillator (5).

6. Proceed with defibrillation/pacing according to defibrillator manufacturer's operating instructions. Make sure that the electrodes are not used for longer (6) than is stated on the packaging.
7. Disconnect electrode connector from defibrillator.
8. After use, carefully remove the electrode using one hand whilst supporting the skin beneath it with the other. Lift the electrode at the end where the electrode cable is attached and slowly peel it off. Ripping off the electrode or peeling it off quickly may damage the skin. Exercise particular caution if the skin is exceptionally sensitive, especially in children but also in elderly patients, diabetics, or patients on a prolonged course of certain medications known to cause dermatological side effects.
9. Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national authority.

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QCS International LTD., Suite 9 Cumberland Business Park, Wardpark Road, Cumbernauld, G67 3JZ Scotland, United Kingdom