EN	NEUTRAL ELECTRODE CABLE
DE	NEUTRALELEKTRODEN-KABEL
FR	CÂBLE POUR ÉLECTRODE NEUTRE
IT	CAVO PER ELETTRODI NEUTRI
BG	КАБЕЛ ЗА НЕУТРАЛЕН ЕЛЕКТРОД
CS	KABEL NEUTRÁLNÍ ELEKTRODY
DA	NEUTRALELEKTRODEKABLER
EL	ΚΑΛΩΔΙΟ ΟΥΔΕΤΕΡΟΥ ΗΛΕΚΤΡΟΔΙΟΥ
ES	CABLE DE ELECTRODOS NEUTROS
ET	NEUTRAALELEKTROODI KAABEL
FI	NEUTRAALIELEKTRODIN JOHTO
HR	KABEL NEUTRALNE ELEKTRODE
HU	NEUTRÁLIS ELEKTRÓDÁKHOZ VALÓ KÁBEL
LT	NEUTRALAUS ELEKTRODO KABELIS
LV	NEITRĀLĀ ELEKTRODA KABELIS
NL	NEUTRALE ELEKTRODEKABEL
PL	KABLE ELEKTROD NEUTRALNYCH
PT_{EU}	CABO DE ELÉTRODO NEUTRO
PT_{BRASIL}	CABO DE PLACAS TERRA
RO	CABLU ELECTROD NEUTRU
SK	KÁBLE NEUTRÁLNYCH ELEKTRÓD
SL	KABEL ZA NEVTRALNE ELEKTRODE
SV	NEUTRALELEKTRODKABEL
KO	중립 전극 케이블
NO	NØYTRALELEKTRODEKABEL
RU	КАБЕЛЬ ДЛЯ НЕЙТРАЛЬНЫХ ЭЛЕКТРОДОВ
TH	สายเคเบิลอิเล็กโทรดที่เป็นกลาง
TR	NÖTRAL ELEKTROT KABLOSU
ZH	中性电极电缆
AR	كابل الإلكترودات المحايدة

Leonhard Lang GmbH, Archenweg 56, 6020 Innsbruck, Austria;
Fax +43 512 392210, medical@LeonhardLang.at

CE

93543 10082975_c 15.05.2024

EN

NEUTRAL ELECTRODE CABLE

Please read this document and keep it safe. Ensure that anyone using this product is familiar with and understands the information contained in this document.

INTENDED USE / INTENDED PURPOSE

Neutral electrode cables are intended to be used as accessories for electrosurgical devices. Neutral electrode cables are non-sterile, reusable, intended to be used only on intact (uninjured) skin and lead energy from the patient to the electrosurgical device.

Indications:

- Neutral electrode cables return high-frequency current from the patient back to an electrosurgical device and are indicated for all monopolar surgical procedures in which high-frequency electrical current is used to cut or coagulate tissue.

Contraindications:

The neutral electrode cable is a medical device that must not be used if one of the following conditions is present:

- Damage to or soiling of the neutral electrode cable that could lead to increased electrical resistance, inadequate contact, or short-circuit.
- Incompatibility of the neutral electrode cable with the high-frequency electrosurgical device or the neutral electrode that could lead to a malfunction or a hazard for the patient or user.
- Use of flammable gases or liquids containing solvents during the surgical procedure in the immediate vicinity of the neutral electrode cable.
- If the surgeon considers that there are additional contraindications or that situations could arise based on the current specialist literature that may endanger the patient or user.

WARNING

- Improper use can cause damage to tissue. These instructions for use are intended for patient safety. **NOT FOLLOWING THESE INSTRUCTIONS MAY LEAD TO BURNS OR OTHER SKIN DAMAGE DURING USE.**
- The neutral electrode cable is non-sterile.
- The neutral electrode cable must not be sterilized.
- Do not use the neutral electrode cable if it is damaged or modified.

PRODUCT LIMITATIONS:

- The neutral electrode cables are only intended for use with Skintact® neutral electrodes.

COMPATIBILITY

The neutral electrode cables are compatible with the dynamic contact quality monitoring systems REM™, NESSY®, and ARM™ and Skintact® neutral electrodes.

CAUTION

- Do not use the neutral electrode cable if the packaging is damaged.
- If the coagulating or cutting effect decrease during surgery or a higher than normal power setting is required, inspect all components of the neutral electrode cable (clamp, cable, connector) as well as the active accessories and the generator before increasing the power setting!
- If you reposition the patient, check the cable connection afterward.
- The neutral electrode cable may be reused. Before using the neutral electrode cable again, it must be disinfected to prevent the risk of a cross-infection from one patient to another.
- Use the lowest effective safe power setting to achieve the desired surgical effect.
- The electrosurgical devices used must be of the type BF (Body Floating) or CF (Cardiac Floating).
- The approved length of accessories defined for the specific electrosurgical device must be compared to the length of the neutral electrode cable. The length of the neutral electrode cable can be found on the packaging information.

USER GROUP

Neutral electrode cables must be applied to the patient by a physician or appropriately trained medical professionals.

INSTRUCTIONS FOR USE

INSPECTION:

Visual inspection of the neutral electrode cable, including the clamp and connector:

- Check for soiling, clean if necessary, and then disinfect (see CLEANING AND DISINFECTION).
- Check for damage (cracks, breaks, discoloration, legibility of the flag label, etc.).

FUNCTION TEST:

- Check the measured signal for CQMS of the neutral electrode cable that is inserted.
- If the measured signal is absent or faulty, check if the connector and clamp are correctly connected.
- If the measured signal is still absent or faulty, replace the neutral electrode cable.

CLEANING AND DISINFECTION:

Cleaning:

- Aqueous sodium hypochlorite solution (10% w/w) is recommended for cleaning.

1. Apply the recommended cleaning solution to a clean cloth.
2. Clean the outside of the neutral electrode cable along the cable insulation.
3. Clean all sides of the clamp of the neutral electrode cable.

Disinfection:

- Only use alcohol-based solutions for disinfecting! Isopropyl alcohol 70% in water and ethanol 70% in water are preferred. Materials may be damaged if the concentrations specified are exceeded!
- The neutral electrode cable can become brittle if allowed to soak in the alcohol-based solution for too long.
- Do not use other disinfection methods (e.g., autoclaving or disinfection baths) because they can damage the neutral electrode cable.
- Do not use solutions containing iodine because this can lead to corrosion of the contact elements.
- Do not sterilize the neutral electrode cable!

1. Apply the recommended disinfectant solution to a clean cloth.
2. Wipe over the surface of the neutral electrode cable, the accessible parts of the opened clamp, and the connector, and ensure that the components are adequately wetted.
3. Close the clamp, wipe over the surface of the closed clamp, and ensure that it is adequately wetted.

Drying:

Allow the neutral electrode cable to air-dry.

APPLYING THE NEUTRAL ELECTRODE CABLE

Connecting the neutral electrode cable:

1. Before connecting the neutral electrode cable, carry out an inspection and a function test (see the sections INSPECTION and FUNCTION TEST).
2. Open the clamp.
3. Insert the contact lug of the adhered neutral electrode completely into the contact opening of the clamp.
4. Close the clamp again by pushing the locking lever down as far as it will go.

5. Check that the connection lug of the neutral electrode is fully inserted into the clamp and does not touch the patient's skin. The clamp must not end up underneath the patient.
6. Position the neutral electrode cable so that the cable does not pull the neutral electrode off the patient's skin. Position the neutral electrode cable to prevent contact with the patient and other cables and to ensure that there are no loops. To prevent capacitive coupling and resulting burns, do not wrap the neutral electrode cable around any of the patient's extremities or other grounded objects.
7. Fix the crocodile clamp at a suitable location to ensure that there is no strain in the cable.
8. Connect the connector to a suitable electrosurgical device.

Disconnecting the neutral electrode cable:
1. Disconnect the connector from the electrosurgical device.
2. Release the crocodile clamp.
3. Open the clamp.
4. Remove the contact lug of the adhered neutral electrode from the clamp.
5. Close the clamp again.
6. Clean and disinfect the cable (see the section CLEANING AND DISINFECTION).

DISPOSING OF THE NEUTRAL ELECTRODE CABLE
- Dispose of the product in accordance with local regulations or the guidelines of the hospital/health care facility.

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Medical Devices Regulation (EU) 2017/745):
If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

Environmental conditions: +10°C to +40°C; 20% to 85% RH

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